



DIGITAL SPECIALTYS

Office Address:-

Shop No-3, Block-D/1, Rajendra Vihar Apartment
Forest Park, Bhubaneswar-751009, Odisha, India

Factory Address:-

Plot No-359/2436, Ground Floor
Near Sukanti Devi Enclave, AIIMS Nagar Lane-1
Patrapada, Aaiginia, Bhubaneswar-751019, Odisha, India

PRODUCT SPECIFICATION

PRODUCT NAME : NITRILE EXAMINATION GLOVES NON STERILE POWDER-FREE

Special Features

- Latex-free, odour-free and powder-free
- Can withstand IPA washing (Iso Propyl Alcohol).
- Can be used in food industry
- Ideal for doctors with latex allergy

Product Information

- Material : Nitrile Rubber, This product does not contain protein found in Natural Rubber Latex
- Colour : Blue
- Design : Ambidextrous, Non-Sterile, beaded cuff and extra rough (textured) at fingertip surface.
- Powder Content : $\leq 2.0 \text{ mg/dm}^2$
- Sterilization : Non-Sterile
- Labeling : Shall comply with the corresponding labeling specification and customer requirement.
- Shelf Life : 3 years from the date of manufacture
- Storage condition : Shall be stored in cool dry place and away from direct light

Packaging Details

Inner Box	100 pcs
Master Carton	2000 pcs (100 x 20 inner boxes)

Certifications : CE
: ISO 9001 : 2015
: ISO 13485 : 2016

Physical Dimensions

Length(mm) Min	Size	XS	S	M	L	XL
	Length	Min 220mm			Min 230mm	
Width (mm)	Size	XS	S	M	L	XL
	Width (+3)	76+3	84+3	94+3	105+3	113+3
	Palm	Min 0.05mm (Single Wall)				
	Finger	Min 0.05mm (Single Wall)				

Office Address:-

Shop No-3, Block-D/1, Rajendra Vihar Apartment
Forest Park, Bhubaneswar-751009, Odisha, India

Factory Address:-

Plot No-359/2436, Ground Floor
Near Sukanti Devi Enclave, AIIMS Nagar Lane-1
Patrapada, Aaiginia, Bhubaneswar-751019, Odisha, India

Physical Properties

Characteristics	Before Ageing	After Ageing (70± 2°C for 166 ± 2 hrs)
Tensile Strength (Mpa) min.	14 mpa Min	14 mpa Min
Ultimate Elongation (%) min.	500% Min	400% Min

Directions for use:

1.1 Preparation For Use : Ensure selection of gloves with correct size.

Below is the correct glove size width:

Size	XS	S	M	L	XL
Width in mm (+3)	76+3	84+3	94+3	105+3	113+3

1.2 How To Use

- Dry hands before donning the gloves.
- Take out the product after checking package.
- Gloves need to be donned carefully; do not apply excessive force to pull the gloves when donning.
- Keep the fingernails short and change gloves frequently to avoid tearing issue.
- The gloves should be discarded with an appropriate waste disposal container after use according to local disposal regulations.

1.3 Cautions For Use

- The product is for single use only.
- Does not contain natural rubber latex.

1.4 Storage Condition

- The atmosphere must always be conducive to the storage of the gloves.
- Gloves before and after opened for use shall be stored in dark, cool and dry place.
- Gloves need to avoid area with excessive heat. The storage areas need to be controlled below 40°C/ 104°F.
- Gloves need to be stored at the areas with good air ventilation.
- Gloves need to be shielded from direct sunlight and fluorescent lighting.
- Gloves need to be stored far away from any energy radiating and heat generating equipment.

1.5 Expiration Date

- The shelf life of the gloves is 3 years from the date of manufacture.